

## Product News

### Advanced Film has Strong Adhesive for PCR Sealing and Storage

**Victorville, CA. February 8, 2008.** New ThermalSeal A™ sealing films are an advanced version of classic ThermalSeal® polypropylene films from Excel Scientific. ThermalSeal A films incorporate a stronger and thicker adhesive layer for more secure sealing of microplates in PCR and storage applications. Although shaped to fit within the rim of raised-rim plates for more reliable sealing of all wells, ThermalSeal A films can also be used for sealing standard microplates. Each film has two removable end tabs that aid in positioning to the plate. Recommended temperature range is -40 °C to +125 °C.

**[www.excelscientific.com](http://www.excelscientific.com)**

"ThermalSeal" is a registered trademark, and "ThermalSeal A" is a trademark of Excel Scientific, Inc.

#### **Contact:**

Excel Scientific, Inc.  
18350 George Blvd.  
Victorville, CA 92394

Phone: (760) 246-4545  
Fax: (760) 246-5588  
Email: [info@excelscientific.com](mailto:info@excelscientific.com)

*An Innovative and Responsive Manufacturer of Laboratory Products*