

Product News

Adhesive Films in Roll Format for Automated Microplate Sealing

Victorville, CA. April 28, 2010. Excel Scientific introduces Roll-Seal™ adhesive sealing films on rolls for use with high-throughput automated microplate sealers. Constructed on three-inch plastic cores, Roll-Seal rolls are compatible with most common adhesive sealers. The Roll-Seal format provides reliable, efficient sealing at a lower cost per plate than sheets or heat-seal films with minimal user intervention. Currently offered are three of Excel's extensive line of adhesive sealing films: ThermalSeal RTS™ clear films with ultra-strong silicone adhesive for qPCR and sitting-drop protein crystallization; AlumaSeal® pierceable aluminum foils for PCR, HTS, and cold storage; and breathable AeraSeal™ films for cell and tissue culture.

www.excelscientific.com

"ThermalSeal" and "AlumaSeal" are registered trademarks. "Roll-Seal", "ThermalSeal RTS", and "AeraSeal" are trademarks of Excel Scientific, Inc

Contact:

Excel Scientific, Inc.
18350 George Blvd.
Victorville, CA 92394

Phone: (760) 246-4545
Fax: (760) 246-5588
Email: info@excelscientific.com

*The Plate-Sealing and Reagent-Handling Experts*SM