

Product News

New Aluminum Foil Film has Special Adhesive for Cold Storage

Wrightwood, CA. November 29, 2006. Excel announces AlumaSeal CS™ adhesive aluminum foil films specially formulated with a new adhesive to withstand long-term storage down to -80 °C. Aluminum foils have excellent barrier properties to retard evaporation. AlumaSeal CS films have a single non-perforated end tab to simplify application. Simply hold the tab and strip the backing from the body of the film as it lays along the plate to avoid curling. Dimensions are 82.6 by 132.6 mm, including the 9.5 mm end tab. AlumaSeal CS films are certified DNase-, RNase-, and nucleic-acid-free. They are available non-sterile or sterile.

www.excelscientific.com

"AlumaSeal CS" is a trademark of Excel Scientific, Inc.

Contact:

Excel Scientific, Inc.
P.O. Box 476
Wrightwood, CA 92397

Phone: (760) 249-6371
Fax: (760) 249-6395
Email: info@excelscientific.com

An Innovative and Responsive Manufacturer of Laboratory Products

P.O. Box 476 • Wrightwood, CA 92397 • Tel: 760.249.6371 • Fax: 760.249.6395 • www.excelscientific.com