


Aspirator Pipets




- for tissue culture media and other aspiration applications
- unique smooth one-piece polystyrene construction
- three sizes and selection of tip styles to suit application
- supplied sterile, individually wrapped in one-side-paper one-side-plastic easy to open packaging

Sterile individually-wrapped aspirator pipets are uncalibrated, with no volume marks and no cotton plugs, providing ease and economy for tissue culture media and other aspiration applications where volume measurement is not required. They eliminate breakage and safety concerns that arise when glass Pasteur pipets are used.



Copyright © 2008. Excel Scientific, Inc.

Catalog No.	Description	Length (inches)	Graduation Interval (mL)	Negative Graduations (mL)	Tip I.D. (mm)	Quantity Per Pkg	Quantity Per Case
ASP-0020-S01	2 mL Standard	10.67	N/A	N/A	1.8-2.2	500	2,000
ASP-010W-S01	10 mL Standard	11.57	N/A	N/A	1.6-2.4	200	800



EXCEL Scientific, Inc.
The Plate-Sealing and Reagent-Handling Experts

18350 George Blvd.
Victorville, CA 92394

email: info@excelscientific.com
760-248-4545 • FAX: 760-246-5588

www.excelscientific.com

"Committed to Exceeding Expectations!"

LINK
www.excelscientific.com/aspirators_content.html