

"Shorties" Serological Pipets

- **much shorter than standard serological pipets for the ultimate convenience when working under hoods**
- **smooth one-piece polystyrene construction eliminates sample hang-up and accidental breakage**
- **superior accuracy: $\pm 1.5\%$**
- **supplied sterile, individually wrapped**
- **one-side-paper one-side-plastic easy-to-open packaging**


Excel's 5 mL and 10 mL "Shorties" serological pipets are significantly shorter than the standard models, for ease of use in really tight spots such as under bench-top hoods. Constructed of virgin polystyrene, they are certified pyrogen-free and non-cytotoxic. Color-coded woven polypropylene plugs provide easy size identification. Shorties are available with standard tips only.

Unique one-piece construction with no seams at either end provides protection against sample hang-up or accidental breakage due to irregularities or weakness. The one-piece construction also provides improved accuracy with calibration to deliver (TD) at $\pm 1.5\%$.

Shorties serological pipets are supplied sterile, individually wrapped in static-free, one-side-paper one-side-plastic packaging that can be opened by peeling the end or by pushing the pipet through the paper side.

Copyright © 2008, Excel Scientific, Inc.

| Catalog No. | Description | Length (inches) | Graduation Interval (mL) | Negative Graduations (mL) | Tip I.D. (mm) | Quantity Per Pkg | Quantity Per Case | |
|----------------|----------------|-----------------|--------------------------|---------------------------|---------------|------------------|-------------------|-----|
| Shorty! | | | | | | | | |
| SER-S050-S01 | 5 mL "Shorty" | | | | | | | |
| | | | 8.66 | 0.10 | 3.0 | 1.6-2.4 | 200 | 800 |
| Shorty! | | | | | | | | |
| SER-S100-S01 | 10 mL "Shorty" | | 8.66 | 0.20 | 8.0 | 2.3-3.5 | 100 | 400 |



EXCEL Scientific, Inc.
 The Plate-Sealing and Reagent-Handling Experts

18350 George Blvd.
 Victorville, CA 92394

email: info@excelscientific.com
 760-248-4545 • FAX: 760-246-5588

www.excelscientific.com

"Committed to Exceeding Expectations!"

LINK

www.excelscientific.com/shorties_content.html